


Trouble hearing?

Do you have trouble hearing in some situations?



At home. Does your partner often need to speak up to get her point across?



With friends. Do you have trouble following the conversation when you are out with friends?



In theaters. Do you avoid going to the movies or concerts?



In stores. Is it hard to understanding what the clerk says when you are out shopping?



Watching TV. Do you miss words when you watch TV at normal volume?



With hearing aid. Do you feel that your hearing aid doesn't fully help in some situations?

This is the solution

The Mino is a new speech amplifier that helps you hear better in most everyday situations, whether you wear a hearing aid or just suffer the occasional hearing difficulty.



Sound with CD quality

The Mino uses similar digital technology as modern CD/DVD players and delivers crystal clear speech, free from background noise.

Pocket size and super light

The Mino is as small as a modern cellphone and weighs just 2½ oz. So it can be kept in your pocket without weighing you down.



Discover Mino



Features.

Packed with useful features.




Single touch action

The Mino is intuitive to use. All functions are accessed with the simple press of a button and there are no complicated menus to navigate.

Long lasting power

The Mino uses a rechargeable lithium-ion battery, just like a modern cell-phone. So you never need to worry about buying new batteries.




Zoom the sound

The Mino has both an omnidirectional and a unidirectional microphone (zoom). With the zoom microphone, you can 'point' to the person talking, reducing any background noise.



Bring Mino to the movies

The Mino is equipped with a telecoil. The telecoil picks up sound in venues fitted with hearing loop systems, such as theaters, churches, and concert halls.



Customize your Mino

If you'd like to tailor your solution, we offer a number of accessories, ranging from car chargers and travel cases to advanced hearing accessories. It's your decision.



Enjoy music and TV

With the Mino, you can enjoy music and TV at any volume without disturbing others. Always with full hi-fi quality.



Boost your hearing aid

Connect a neck loop, and the Mino further amplifies speech and acts as an additional volume and tone control for your hearing aid.



MicroSet™

A discreet way to better hearing.

How it works.



Discover how the Mino can help you in your daily life.

At home



Listen at the volume of your choice without others having to raise their voices. Simply place the Mino near the other person and worry no more.

The Mino effectively isolates speech while removing disturbing background sounds. The result is crystal clear speech with no noise. We offer a number of accessories for the Mino. Non-hearing aid users can, for example, choose to listen using headphones, MicroSet™ earphones or steto-clips.



With friends



Trouble hearing people in noisy environments? Change to the zoom microphone to reduce background noise and focus on the conversation.

The zoom microphone is an excellent choice in environments that features a lot of background noise like for instance restaurants, coffee houses, train terminals, stores or airports. You can easily change between the microphones with a simple press of a button, depending on your current need.

In theaters



Do you avoid movie theaters and concerts? Many public venues, such as theaters, concert halls, and churches, are fitted with hearing loops.

Sound from the loudspeakers in such venues is converted into an electromagnetic field that the Mino picks up via its built-in Telecoil. The Mino's Telecoil always delivers perfect sound, whether it's placed in your pocket or on a table. Look for a T symbol sign next time you visit your favorite movie theater.



Multitalented.



Discover how the Mino can help you in your daily life.

In stores



When you're shopping down town, it's handy to keep the Mino in a pocket so that both your hands are free.

The MicroSet™ is a unique accessory specially developed for such times. It includes a small lapel microphone that you simply attach to your clothes and a discreet earphone that you pop in your ear. The lapel microphone is unidirectional (zoom) and reduces background noise in the store, so you can talk with the clerk uninhibited.



Watching TV



When you have difficulty hearing, it's easy to turn the TV volume up to the max, which other family members don't always appreciate.

With the Mino, there's no reason to turn up the TV volume to hear better. Use a cable kit to connect the Mino to the TV and you can enjoy great TV sound in stereo while sparing your family's ears. If you want to comment on the show, simply press the Push-to-Listen button on the Mino to lower the volume and join the conversation.


With hearing aid



If your hearing aid isn't quite enough in certain situations, you should try connecting it to your Mino using our discreet neck loop.

The Mino further amplifies the sound and also acts as a remote for easy adjustment of the volume and treble, as well as microphone selection. You can also connect the Mino to the TV to amplify the TV sound in your hearing aid.

Do you travel a lot? Then you can store your Mino and accessories in this protective travel case for the Mino.



Accessories.

Create a solution completely tailored to your needs.



Headphones BE9122
Robust stereo headphones with great sound quality. Features semi closed drivers that reduce sound leakage. Foldable for easy travel.



Stetoclips BE9125
The Stetoclips are capable of delivering sound pressure up to 140 dB. They are comfortable to wear and very easy to put on and take off.



Earphones BE9124
Small and light weight stereo in-ear earphones with great sensitivity and excellent isolation from ambient noise.



Neckloop BE9128
With a neck loop, the Mino becomes best friend with your hearing aid. Just plug it into your Mino and switch your hearing aid(s) to T-position.



Cable kit BE9137
Use the cable kit to connect the Mino to your TV or multimedia equipment and enjoy high quality stereo sound. Tele plug and RCA adapter is included.



Tie-clip microphone BE9136
Attach it close to the sound source in difficult listening environments. Also a great choice for a discreet look as the Mino can be hidden in a pocket.



USB car charger BE9143
With the car charger you can easily charge your Mino unit in your car. Simply plug it in, and your Mino will be fully charged within 2.5 hours.



Travel case BE9147
The travel case is made of black and white canvas with reinforced cover. It effectively protects the Mino and your accessories during transport.



MicroSet™
The MicroSet includes a small directional lapel microphone that you can attach to your clothes and a discreet earphone that you pop in your ear.



Bellman & Symfon®

Head Office
Bellman & Symfon Europe AB
Södra Långebergsgatan 30
421 32 Västra Frölunda
Sweden
Phone +46 31 68 28 20
E-mail info@bellman.se
www.bellman.se

China Office
Bellman & Symfon Asia Ltd.
1115-16, Main Tower, GD Int'l Hotel
339 Huanshi Dong Rd
Guangzhou, China, 510095
Phone +86 20 8318 2950
E-mail info.cn@bellman.com
www.bellman.cn